

Cole Relaxation Frequency: A Parameter for Pancreatic Cancer Detection

¹Les Bogdanowicz,¹Martina Guidetti, ¹Onur Fidaner, ¹Daniel Gehrke, ¹Alexander Gryczuk, ¹Donato Ceres, ¹Margaret John, ¹Constantine Bovalis, ¹Erik Kundro, ²Karla Castellanos, ²Akshar Patel, ³Isaac Rajiman, ²Paul Grippo

¹NovaScan Inc, Chicago, IL, USA; ²University of Illinois at Chicago, Chicago, IL, USA, ³Texas International Endoscopy Center, Houston, TX, USA

INTRODUCTION

- While there are multiple tools available to aid in Pancreatic Cancer (PC) diagnosis including imaging methods such as CT, MRI, ERCP and endoscopic ultrasound, real-time, on-the-spot diagnosis can be challenging.^{1,2}
- Lesions are generally sampled via biopsy, which presents several disadvantages; lesion heterogeneity, inter-observer variability in identifying and grading the lesions, low diagnostic yield due to insufficient integrity or size of samples, complications, and inconclusiveness of the initial biopsies.^{1,2}
- There is demand for an on-the-spot, real-time assessment device that works as a **decision support tool** for the endoscopist without image interpretation.

OBJECTIVE

- To determine if Novascan technology based on Cole-Relaxation Frequency (CRF) calculation from impedance spectroscopy can detect cancer in pancreatic tissues.
- To demonstrate if CRF can predict the level of fibrosis in PC


METHODS

- The genetically engineered mouse model (GEMM): *LSL-Kras;LSL-p53;Pdx1-Cre* (KPC) was used as a platform for investigating CRF in PC
- CRF was measured *ex vivo* with NovaScan technology on excised pancreas samples in a cohort of n=26 genotyped mice, of which
 - 15 were KPC
 - 2 were KC (model of pancreatic neoplasia)
 - 9 were wild type controls

METHODS cont'd

- Cancer presence is determined if **CRF > 1 MHz**
- Outcomes were compared to histopathology results for each of these samples to determine the specificity and sensitivity of the NovaScan device
- Pancreases were scored as percent fibrosis over multiple fields of view and compared to CRF
- To determine if CRF can effectively discern fibrosis in PC from pancreatitis, n=18 cerulean-induced acute pancreatitis (AP) and n=6 saline injection controls were also tested


RESULTS




- NovaScan technology identified 4 out of 5 PanIN samples as cancerous. **Specificity** and **sensitivity** were **100%** and **94%**, respectively on the entire cohort and **100%/100%** after excluding the PanIN samples.
- Cerulean-induced AP and Saline injected controls tested negative.

CONCLUSION

- NovaScan's Impedance Spectroscopy CRF technology discerns healthy and cancerous lesions in pancreas, similarly as shown in breast skin and lung previously.^{3,4,5}
- The CRF can predict the level of fibrosis in PC, results are similar to those shown in breast cancer previously.³
- As the CRF value increases the likelihood of cancer increases.



- Development of this technology includes a digital rapid-on-site evaluation, "**digital-ROSE**" for **quick assessment** of biopsy samples



References:

- cancerstatisticscenter.cancer.org
- Zhang et al., Challenges in diagnosis of pancreatic cancer, 2018.
- Gregory et al., The Cole relaxation frequency as a parameter to identify cancer in breast tissue, 2012.
- Svoboda et al., Bioimpedance measurement as an assessment of margin positivity in Mohs surgical specimens of nonmelanoma skin cancer: Management implications, 2018
- Bogdanowicz et al., The Cole Relaxation Frequency as a Parameter to Identify Cancer in Lung Tissue: Preliminary Animal and Ex Vivo Patient Studies, 2022..